Manuel d'instructions
Bedienungsanleitung
Manuale di istruzioni
Manual de instrucciones

PERKINS MK3 TONOMETER

Instructions for use

For full instructions for use and information on all products in the Haag-Streit UK portfolio, please visit **www.haagstreituk.com**

Contents

en	4
fr	34
de	64
it	94
es	124

Introduction

Read the instructions for use carefully before commencing use of this device. It contains important information regarding the safety of the user and patient.

Applanation tonometry may only be carried out by appropriately trained and qualified medical personnel.

Device Description

The Perkins tonometer is an ophthalmic, battery-powered, reusable measuring instrument with reusable test weights and either a reusable prism or a tray of 20 Tonosafe disposable prisms.

Intended Use

The Perkins MK3 tonometer is intended to measure Intraocular Pressure (IOP) by applanation (applying a small flat disc to the cornea).

The Perkins Tonometer Test weights are used to apply a predetermined force to the Perkins tonometer to show confirmation that the readable gauge displays the same value as the weight.

Intended purpose

Perkins MK3 Tonometer:

An ophthalmic, battery-powered, reusable measuring instrument with reusable Test weights to be used with a single use or limited reuse (100 processing cycles or 2 years) prism. It is intended for use by healthcare professionals for determining the intraocular pressure (IOP) by exerting an external force on the cornea to which the resistance deformation is measured in millimetres of mercury (mmHg). The Tonometer body, Test weights and prism can contact the user's skin. The Tonometer forehead rest contacts the patient's body, while the prism contacts the patient's cornea of each eye up to three times during the procedure. The total patient exposure is determined as having a cumulative transient patient contact, i.e. less than one minute. The Tonometer and Test weights touch the user with limited cumulative contact. Protective medical gloves are recommended when using the Tonometer giving limited exposure to the user. The Tonometer is able to be used on any patient population and is only contraindicated where medical assessment by the healthcare professional has deemed it unsafe in the case of infection or trauma. The Test weights are only to be used with the Perkins Tonometer.

Perkins Tonometer Test weights:

A reusable medical device accessory intended to be used to provide a known force to the reading scale on the Perkins MK3 Tonometer. Contact is transient use with the user's hands.

Contraindications

The Perkins tonometer is contraindicated where medical assessment by the healthcare professional has deemed it unsafe in the case of infection, disease or trauma.

Medical assessment is necessary and due care is to be exercised in cases of eye infections and injuries to the cornea.

A list of known contraindications is as follows;

- Trauma/corneal ulcer
- Ruptured globe
- Infection
- Unhealed corneal abrasions or ulcers
- Inability to administer topical anaesthetic

Manufacture

The Perkins tonometer is manufactured in the UK by Haag-Streit UK (HS-UK). Information on the Haag-Streit UK product portfolio can be found on the HS-UK website at www.haagstreituk.com

en

Contents

1. Safety	7
2. General description	11
3. EMC declaration	13
4. Appliance assembly	14
5. Operation	16
6. Technical data	25
7. Maintenance	26
8. End-of-life & disposal	30
Appendix A: Accessories, consumables, spare parts	31
Appendix B: Legal	32

1 Safety

Failure to comply with these instructions for use may result in damage to the device and/or pose a danger to patients or users.

The device must not be stacked or placed in close proximity to other electronic devices.

Some portable radio frequency equipment, like mobile phones or RF telephone equipment, including antennas, may cause interference with medical devices. Such equipment has to be kept at a distance of more than 30cm (12 inches) from any part of the instrument. Failure to observe this precaution may result in the device malfunctioning. If the distance cannot be observed, ensure the tonometer is functioning as expected before use.

Only approved Haag-streit UK components should be used with this device and repairs should be carried out by HS-UK trained engineers.

Personal Protective Equipment (PPE) must be worn when using this device

Any serious incident in relation to this device must be reported to the manufacturer and the regulatory authority of the country where the event occurred.

1.1 Areas of application of the device

The device is intended for use in professional healthcare facility environments, including, but not limited to; doctor's practices, hospitals, optometrists and opticians premises.

1.2 Environmental conditions

Transport	Temperature	-40°C to +70°C
	Atmospheric pressure	500hPa to 1060hPa
	Relative humidity	10% to 95%
Storage	Temperature	-10°C to +55°C
	Atmospheric pressure	1060hPa
	Relative humidity	10% to 95%
Use	Temperature	+10°C to +35°C
	Atmospheric pressure	800hPa to 1060hPa
	Relative humidity	30% to 90%

Do not use the device if it is exposed to environmental conditions outside of those specified above. Contact Haag-Streit UK Ltd for further guidance.

1.3 Shipment & unpacking

Do not use the device if it is unintentionally opened before use. In this case, contact Haag-Streit UK Ltd for further guidance

Check the appliance for damage after it is unpacked. Do not use damaged equipment. Return defective appliances in the appropriate packaging to HS-UK for repair.

Retain original packaging. Use original packaging when transporting or returning for servicing.

The device should be cleaned after unpacking before its first use.

1.4 Installation warnings

Only original manufacturer-supplied replacement parts and materials may be used.

Do not modify this equipment.

Repairs may only be performed by fully-trained Haag-Streit UK accredited engineers. Contact the Haag-Streit UK Service Division for a quotation. Contact details are available at www.haagstreituk.com/service. Users who live outside the UK should contact their local Haag-Streit UK representative.

1.5 Operation environment

This equipment must only be operated by qualified medical professionals. The clinical institute is responsible for any training on the device. The Perkins tonometer may only be used in accordance with the intended use.

Never use the device in a potentially explosive environment where volatile solvents (alcohol, petrol etc) and flammable anaesthetics are in use.

The examination is to be performed at the medical professional's discretion and consideration should be taken in cases of active infections and lesions of the eye. Certain conditions could falsify measured values and the procedure may interfere with the patient's condition.

Strong magnetic fields may falsify the measurement results.

Use only original Haag-Streit Goldmann tonometer measuring prisms or the Tonosafe sterile disposable prisms available from Haag-Streit UK.

It is known that refractive surgery has an effect on Goldmann Applanation Tonometry (GAT) values, as the procedure changes the biomechanical properties of the cornea. This could cause an underestimation of the pressure and therefore, secondary checks are required to confirm diagnosis.

The biomechanical properties of the cornea have an effect on the measurement readings.

1.6 Tonosafe

Tonosafe was developed to avoid the necessary cleaning and disinfection of measuring prisms from applanation tonometers. Tonosafe is a disposable applanation and optical duplication prism for Goldmann and Perkins applanation tonometers, for single use.

Additional details can be found in the separate Tonosafe instructions for use.

1.7 Disinfection

Do not re-use Tonosafe sterile disposable prisms.

The Haag-Streit Goldmann tonometer measuring prisms are made of Polymethyl methacrylate (PMMA). Please take note of the following restrictions:

- Do not disinfect with alcohol
- Do not clean with acetone
- Do not disinfect using UV radiation
- Do not sterilise using steam or ethylene oxide
- Do not expose to temperatures above 60°C

The Haag-Streit Goldmann measuring prisms are not shipped disinfected and must be cleaned and disinfected before their first use in accordance with the separate instruction manual on the cleaning and disinfection of tonometer measuring prisms and contact glasses.

- The operator shall be liable in the event of non-observance of the cleaning and disinfecting process.

The above-mentioned instruction manual, the brief instructions, information on the disinfection of tonometer measuring prisms and contact glasses and additional information can always be found on the website www.haagstreituk.com


1.8 Warranty & product liability

Haag-Streit UK products must be used only for the purposes and in the manner described in the documents distributed with the product.

Improper handling can damage the product. This would void all warranty claims.

Continued use of a damaged device may lead to injury. In such a case, the manufacturer will not accept any liability.

1.9 Description of symbols

	Name, address and contact details of the manufacturer		Great Britain Conformity Assessment for medical devices
	Swiss authorised representative		Consult instructions for use
	Product Identifier		Caution, consult accompanying documents
	Country and date of manufacture		Read instructions manual
	LOT number		Type B applied parts
	Serial number		Temperature limit
	Unique device identifier		Humidity limitation
	Quantity		Atmospheric pressure limitation
	EU Authorised Representative		Do not dispose of with household refuse
	Medical device		Do not use if packaging is damaged
	Conformité Européene, European conformity assessment for medical devices		Keep Dry
	Direct Current		Recycle
	Keep out of sunlight		

2 General Description

2.1 Overview


Figure 1

1. Forehead rest
2. Forehead rest locking screw
3. Forehead rest stem
4. Prism holder
5. Doubling prism
6. Illumination LEDs

Figure 2

- 7. Viewing lens
- 8. Milled thumb-wheel
- 9. Scale
- 10. Tonogrip (Cone arm securing buttons)
- 11. Battery Indicator
- 12. Charging Port

3 Electro Magnetic Compatibility (EMC)

When in use, medical devices may be susceptible to electronic interference from other electronic instruments, such as PCs and mobiles. This may result in incorrect operation of the medical device and create a potentially unsafe situation and possible EMC shock. Medical devices should not interfere with any other devices.

In order to regulate the requirements for EMC with the aim to prevent unsafe product situations, the EN 60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices.

Nevertheless, special precautions need to be observed:

The use of accessories and cables other than those specified by Haag-Streit-UK Ltd, with the exception of cables sold by Haag-Streit-UK Ltd as replacement parts for internal components, may result in increased emission or decreased immunity of the device. Only use HS-UK replacement parts for this device.

The medical device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is unavoidable, the medical device should be checked prior to use for normal operation.

Further guidance regarding the EMC environment (in accordance with EN 60601-1-2) in which the device should be used is available at:

www.haagstreituk.com

4 Appliance installation

4.1 Set up

Remove the Perkins tonometer from the protective case by holding the main body. Under no circumstances should it be lifted by the prism holder.

Figure 3

1. Perkins Tonometer Mk3
2. Forehead Rest
3. Doubling Prism/Black Tonosafe Holder
4. Two-gram weight

5. Five-gram weight
6. Battery Charger
7. Adaptor (UK, US, EU)
8. Setting Block

5 Operation

Full operational check of all functions must be carried out prior to use. If there is any doubt, do not use the device and contact Haag-Streit UK.

Differences in atmospheric pressure can have an effect on results and therefore, secondary checks should be carried out to confirm diagnosis.

Fluorescein sodium is required in this procedure, ensure user has no allergies.

Regular calibration is required to maintain correct functioning of the device. Refer to Section 7 of this IFU.

Do not use without a Tonosafe prism or Haag-Streit Goldmann prism.

Risk of photo-toxicity, restrict exposure to LEDs as much as possible.

Blue LEDs may stop working upon Electrostatic Discharge (ESD) contact discharge to the thumb wheel at +8kV and -8kV. In this instance, the device is not damaged and requires resetting by turning the device off and back on.

Do not put the Perkins Tonometer on a hot surface, damage will occur.

Take care not to apply too much pressure to the eye or eyelid, risk of injury to patient.

If the device is used on animals, the conformity assessment mark is invalidated as the approval is for medical devices for use on human beings.

The reading dial must be reset between uses.

Avoid long periods of use as the battery and charger could become hot.

Consumption of fluids may increase IOP readings.

Be aware of diurnal variations whilst taking readings.

5.1 How pressure is measured

The applanation tonometer functions according to the 'Goldmann method'; measurement of the pressure required to maintain a uniform applanation of the surface of the cornea.

The readings obtained are based on the assumption that a cornea with a 'normal' thickness is being measured (a different corneal thickness leads to changes in the measured IOP). A 'normal' corneal thickness is considered to fall within the range of 530 to 560 Micrometers.

If wrong measuring results are suspected, please conduct a functional verification in accordance with the 'Maintenance' section.

5.2 Prepare the prism for use

The prism tip should be prepared for use by following the manufacturer's instructions

Tonosafe prism tips and the prism holder have the same mass as the Haag-Streit Goldmann prism. Conversion of results is not required.

To insert the prism into the prism holder;

1. Lean the tonometer back until the prism holder moves forward
2. Gently pinch together the Tonogrip (cone arm securing buttons) on either side to hold the prism holder in position
3. With the cone arm buttons pressed and the prism holder out, the prism may be pushed into place.

5.3 Preparing the patient

Ensure the patient is relaxed before starting the procedure. They should have a normal heart rate and should be in the correct (e.g. seated, supine) position >5 mins before taking the reading.

Both eyes must always be anaesthetised (e.g. 2-3 drops each of an anaesthetic within half a minute) to reduce movements of the lids during examination.

Place a fluorescein paper strip near the external canthus in the lower conjunctiva sac. After a few seconds the lacrimal fluid is sufficiently coloured and the paper can be removed.

Alternatively, combined anaesthetic fluorescein eye drops can be used.

When using drops, a solution of sodium fluorescein 0.25% to 0.5% is advisable.

5.4 Instruction to the patient

Make sure that the measuring arm of the tonometer can move freely and is not hindered by the patient's eyelashes, cheek, etc.

Figure 4

Instruct the patient to look straight ahead or slightly upward and if necessary use a fixation light for this purpose.

It is recommended that patients are repeatedly asked to keep their eyes wide open during the examination. If required, the lids of the examined eye may be held open by the examiner's fingers provided that no pressure is applied to the eye.

5.5 Astigmatism

If the cornea is spherical, measurements can be made on any meridian, but it is most convenient to do it on the 0° meridian.

This is not so when eyes with higher corneal astigmatism than 3 dioptres are examined, as the flattened areas are not circular but elliptic.

It has been calculated that, in cases of larger corneal astigmatisms, a surface of 7.354 mm² (Ø 3.06 mm) is to be applanated, when the measuring prism is at an angle of 43° to the meridian of the greatest radius.

For example: If the corneal astigmatism amounts to

$$6.5\text{mm}/30^\circ = 52.0\text{D}/30^\circ \text{ and } 8.5\text{mm}/120^\circ = 40.0\text{D}/120^\circ$$

The graduation value 120° of the prism is set at the 43° mark of the prism holder.

If there is a corneal astigmatism of $8.5\text{mm}/30^\circ = 40.0\text{D}/30^\circ$ and $6.5\text{mm}/120^\circ = 52.0\text{D}/120^\circ$

The graduation value 30° is set at the 43° mark.

In other words, set the axial position of the greatest radius, that is the axis of a minus cylinder, on the prism graduation mark on the prism holder.

5.6 Measurement

The Perkins tonometer should be held so that the thumb rests on the milled wheel (Fig 2, 8) controlling the spring. The light is switched on by turning the thumb-wheel until the scale reading (Fig 2, 9) is just above zero.

If the forehead rest (Fig 1, 1) is to be used, the stem (Fig 1, 3) should be extended after loosening the locking screw (Fig 1, 2). When the forehead rest is correctly positioned, tighten the locking screw and the instrument can then be brought towards the patient so that the face of the prism comes into contact with the cornea.

For ease of use, hold the tonometer obliquely with the handle slanted away from the nose (see section 5.4). Care should be taken to prevent the prism touching the eyelid margin.


When seen through the viewing lens (Fig 2, 7) the semicircles of fluorescein should now be clearly visible and appear equal in size, denoting that the contact point is central on the cornea.

The light is switched off by turning the thumb-wheel until the scale reading (Fig 2, 9) is below zero.

Measurement of the force necessary to flatten a corneal surface of constant \varnothing 3.06 mm size.

Figure 5

Dial indicator: 1
 Force: 9.81 mN
 Pressure: 1.33 kPa = 10mm Hg

Adjust the force by turning the thumb-wheel (Fig 2,8) until the inner-margins of the semicircles coincide (Fig 5).

The tonometer is removed from the eye and the reading noted. The large divisions of the scale represent grams and the small divisions 0.2 grams.

Adjust the readings taking into account the age and biomechanics of the patient's eye.

The reading is multiplied by ten to give the tension in millimetres of mercury (mm Hg). Readings should be repeated until a steady value is obtained; this will usually be slightly lower than the initial reading.

If the semicircles appear large and are not reduced by diminishing the force of the spring, the tonometer has been pushed too close to the eye and withdrawing it slightly will bring the prism within the range of free movement.

If it is desired to use the Perkins tonometer without any forehead rest, the instrument can be steadied by resting the fingers against the patient's malar region while the prism is in contact with the cornea.

5.7 Operating the Instrument

1. Immediately before taking the measurements, the patient should be made to close the eyes briefly so that the cornea becomes sufficiently moistened with the lacrimal fluid and fluorescein.
2. By moving the tonometer forward, the measuring prism is brought into contact with the centre of the cornea over the area above the pupils. The limbus of the cornea will be illuminated with a bluish hue. This illumination is best observed by direct sight from the opposite side of the illumination unit by the examiner.

As soon as the limbus of the cornea illuminates, cease any further forward movement of the tonometer immediately.

3. After obtaining contact, commence observing the cornea. A regular pulsation of the two semi circular fluorescein rings, which can be of varying size dependent on the ocular pressure, when the measuring drum is at setting 1, will show that the tonometer is in the correct measuring position
4. The pressure on the eye is increased by turning the milled thumb wheel on the tonometer, until the edges of both fluorescein rings just meet (Fig. 7). The edges cross over each other with the pulsation of the eye.

The width of the fluorescein ring around the contact position of the measuring prism should be $\frac{1}{10}$ approx. of the diameter of the applanation surface (0.3mm).

5. Reading the scale:

- The reading
- Multiplied by 10
- Gives the ocular pressure in mm Hg

Note

The following are schematically simplified illustrations.


Figure 6

Semi circular images in the centre of the field of sight






Figure 7






Correct final position

5.71 Sources of Error






Flourescein ring / Distance to the patient

Result	Image in ocular	Cause	Correction
The fluorescein ring is too wide	 A circular diagram representing the ocular view. A horizontal line bisects the circle. A thick, dark, semi-circular ring is drawn across the center, extending slightly beyond the horizontal line on both sides.	The measuring prism was not dried after cleaning, or the eyelids came into contact with the measuring prism whilst measuring.	The tonometer should be withdrawn and the measuring prism dried with a cotton wool swab
The fluorescein ring is too narrow	 A circular diagram representing the ocular view. A horizontal line bisects the circle. A thin, light-colored, semi-circular ring is drawn across the center, extending slightly beyond the horizontal line on both sides.	The lacrimal fluid has dried.	Let the patient close the eyes once or twice. Then repeat the measuring procedure.
No semi circular rings appear, only the centre line	 A circular diagram representing the ocular view. A horizontal line bisects the circle. No other features are present.	The measuring prism is not touching the cornea. Should the patient withdraw their head slightly, the irregular pulsations will occur and the prism will only contact the eye intermittently. Should the patient withdraw the head further, then the fluorescein rings will disappear together.	Steady patient's head
Both of the too large semi circular rings appear partly	 A circular diagram representing the ocular view. A horizontal line bisects the circle. Two thick, dark, semi-circular rings are drawn, one above and one below the horizontal line, each extending slightly beyond the horizontal line on both sides.	Should the tonometer be moved forward towards the patient, or should the patient move towards the tonometer whilst the measurement is being taken, then the feeler arm will come into contact with a sprung stop piece. The applanation surface is then too large.	The image will not change when turning the milled thumb wheel. Withdraw the tonometer until the regular pulsations of a corresponding smaller applanation surface indicate the correct measuring position and pressure changes lead to immediate applanation surface changes.





5.72 Inaccurate Horizontal Centering

Result	Image in ocular	Cause	Correction
<p>The upper semi circular rings appear partly</p>		<p>The measuring prism is not focused on the eye. The eye is too far on the right</p>	<p>Move the tonometer to the right</p>
<p>The upper semi circular ring appears completely - the lower ring partly</p>		<p>The measuring prism is not focused on the eye. The eye is still too far on the right.</p>	<p>Move the tonometer to the right</p>
<p>Correct position: The edges of both semi circular rings meet exactly in the centre of the ocular</p>		<p>The measuring prism is focused precisely on the eye.</p>	
<p>The lower semi circular ring appears completely - the upper ring partly</p>		<p>The measuring prism is not focused on the eye. The eye is still too far on the left.</p>	<p>Move the tonometer to the left.</p>
<p>The lower semi circular ring appears partly</p>		<p>The measuring prism is not focused on the eye. The eye is too far on the left.</p>	<p>Move the tonometer to the left.</p>

5.73 Inaccurate Vertical Centering

Result	Image in ocular	Cause	Correction
A semi circular ring appears partly in the upper half		The measuring prism is not focused on the eye. The eye is too far up.	Move the tonometer upwards.
The ring appears completely in the upper half		The measuring prism is not focused on the eye. The eye is still too far up.	Move the tonometer upwards.
The ring appears almost completely in the upper half, and partly cut in the lower half		The measuring prism is not focused on the eye. The eye is still too far up.	Move the tonometer upwards.
Two partly cut rings appear, the larger one in the upper half		The measuring prism is not focused on the eye. The eye is still too far up.	Move the tonometer upwards.
Correct position: The edges of both semi circular rings meet exactly in the centre of the ocular		The measuring prism is focused precisely on the eye.	

5.74 Inaccurate Pressure

Result	Image in ocular	Cause	Correction
<p>The outer edges of the fluorescein rings touch each other</p>		<p>Pressure is too strongly reduced</p>	<p>Increase the pressure by turning the milled thumb wheel</p>
<p>The fluorescein rings coincide and form a line</p>		<p>The pressure is reduced</p>	<p>Increase the pressure by turning the milled thumb wheel</p>
<p>The fluorescein rings do not touch each other</p>		<p>The pressure is too strongly increased</p>	<p>Reduce the pressure by turning the milled thumb wheel</p>
<p>Correct position: The inner edges of the fluorescein rings touch each other</p>		<p>Accurate pressure - precise focusing of the measuring prism</p>	

6 Technical

Measuring range:	0 - 50 mm Hg
Measurement uncertainty of the force impingement on the measuring prism over a measuring range, from	0 - 49.03 mN
Standard divergence:	$0.49\text{mN} \leq 3s \leq 1.5\%$ of rated value
Net weight, packed:	1.0Kg

7 Maintenance

Do not modify this equipment without authorisation of the manufacturer. Repairs may only be performed by fully-trained Haag-Streit UK accredited engineers. Contact the Haag-Streit UK Service Division for a quotation. Contact details are available at www.haagstreituk.com/service. Users who live outside the UK should contact their local Haag-Streit UK representative.

Only original manufacturer-supplied replacement parts and materials may be used.

Servicing & maintenance

For information on the servicing or repair of the Perkins tonometer, users will need to contact the HS-UK Service Division, at service.uk@haag-streit.com. Information on the HS-UK Service Division can be found on the HS-UK website at www.haagstreituk.com/service

7.1 User Calibration check

Calibration of the Tonometer is advised before every session using the following method;

1. Insert the prism tip (Fig 1, 5) into the prism holder

Check for zero position

2. Adjust the milled thumb-wheel (Fig 2, 8) so that the scale reading is below zero by the full thickness of one scale line. Whatever the angle of the instrument in the hand the prism tip should tend towards its backward position.
3. Set the scale reading above zero by the full thickness of one scale line. The prism tip should tend towards its forward position regardless of the angle of the tonometer.

Check at the 2g position

4. To make a rough check adjust the milled thumb-wheel (Fig 2, 8) so that the scale reading is below the 2 mark by the full thickness of one scale line. Place the instrument on a flat horizontal surface with the setting block under the body and the prism upwards as shown in Fig 8. Place the 2g weight, recess downward, centrally on the prism cone and check that it carries the prism positively down to it's lowest position.
5. Also, as a rough check, set the scale reading above the 2 mark by the full thickness of the scale line and place the instrument on the setting block as before. Place the 2g weight on the prism and check that it remains in its highest position.

If desired, the Actual point of balance can be established by adjusting the milled thumb wheel (Fig 2, 8) until the prism floats at the calibration-position (rim of prism holder level with, or not more than 1 mm above the case surface).

Check at the 5g position

- Repeat the process as described in points 4 and 5 above, but using the 5 mark on the scale and a setting one-and-a-half lines width above and then below the 5 mark and with the 5g weight on the prism cone. The actual point of balance can again be found if required.

The tolerances in the tests above represent the following figures

Zero	$\pm 0.05\text{g}$
2 g	$\pm 0.05\text{g}$
5 g	$\pm 0.075\text{g}$

If the tonometer is found to be outside the stated limits it should be returned to HS-UK for repair.

7.2 Repair

To guarantee a long service life, the device must be cleaned as described and protected in its case when not in use or during transportation. HS-UK recommend annual servicing, by a fully-trained Haag-Streit UK accredited Engineer.

Repairs may only be performed by fully-trained Haag-Streit UK accredited engineers. Contact the Haag-Streit UK Service Division for a quotation. Contact details are available at www.haagstreituk.com/service. Users who live outside the UK should contact their local Haag-Streit UK representative.

7.3 Cleaning & disinfection

Turn the device off before proceeding with cleaning and disinfection.

This guidance does not apply to tonometer measuring prisms.

Please follow national and local infection control laws, standards or guidelines.

Please observe the manufacturer's instructions for all cleaning and disinfecting agents.

Haag-Streit UK only approves the following cleaning and disinfecting procedure. Other procedures and cleaning agents must be validated by the user, any damage to the device is not covered under warranty.

Some parts of the device are not designed for frequent cleaning and disinfection, and may become damaged after some time. These include:

- Eyepiece
- Device labels

The device labels are not designed for frequent cleaning and disinfection. These may degrade if the process is repeated multiple times. If labelling of the device has deteriorated the device must be decommissioned immediately

The optics (Fig 2, 7) are not meant for frequent disinfection and may be permanently damaged if the process is repeated many times.

Do not submerge the device and prevent any ingress of liquid, this will damage the device.

Parts of the device that are in contact with the Patient and/or Operator must be cleaned and disinfected using CE/UKCA marked surface wipes, >70% isopropyl alcohol or ethanol alcohol, prior to every examination.

During use, the external optical surfaces of the tonometer may become soiled, particularly the outside surface of the eyepiece lens. Normally this is the only optical surface that may require regular cleaning.

If required, wipe the optics with lens cleaning paper and cleaning solution starting in the centre of the lens, wipe in a circular motion slowly and radially toward the edge of the lens. Allow the lens to dry completely before use.

Safely dispose of cleaning materials in the local designated waste.

7.4 Charging the battery

The battery pack should only be recharged using the charger supplied with the unit.

If the LED alternately flashes between orange and green, this indicates a battery fault. In this case return the tonometer and battery charger to the manufacturer.

Replacement of the battery pack should only be carried out by the HS-UK engineers.

7.4.1 Charger unit

The Perkins Mk3 tonometer is powered by a 3.6V Ni MH chargeable battery pack that will last for approximately 25 hours.

The battery indicator LED indicates;

Green		Fully charged.
Orange		Several hours of use left.
Red		Battery needs to be charged

When the tonometer is first plugged in, the battery charger LED will show orange for several seconds while it gives a short charge to the tonometer battery.


The battery charger will then switch to yellow for several seconds while it carries out initialisation analysis of the tonometer battery. If the tonometer battery condition is acceptable the battery charger LED will switch to orange again for the "fast" charge and stay orange until Tonometer battery is almost fully-charged (up to 3 hours).

When the battery charger detects the tonometer battery is almost fully-charged, the battery charger LED switches to green with short flashes of yellow (up to 1 hour) to complete the charge.

When the top-off charge is complete, the battery charger LED will switch to solid green (trickle charge) indicating the tonometer battery is fully charged

The battery pack requires no maintenance and is designed to last for the life-time of the device (5 years).

8 End-of-life & disposal

Electrical and electronic devices must be disposed of separately from household and clinical waste.

The battery may contain toxic or hazardous material and must be disposed of in the appropriate electronic waste stream (E.G. WEEE Directive). This will guarantee that no hazardous substances will enter the environment and that valuable raw materials are recycled.

For correct disposal of this electrical device, please contact Haag-Streit customer services on +44 (0)1279 414969, who will instruct you how to return the device for disposal.

Appendix

Appendix A1. Accessories/consumables/spare parts

Do not modify this equipment without authorisation of HS-UK. Installation and repairs may only be performed by HS-UK. Contact your Haag-Streit representative for installation, repairs and modification work on the system. The contact details are available at www.haagstreituk.com

Only original Haag-Streit replacement parts may be used;

- 2701411 Carrying case
- 2701084 Five-gram weight
- 2701086 Two-gram weight
- 2701400 Setting block
- 2701417 Patient's forehead rest
- 2700701 User's Instruction Manual
- 7220316 Doubling prism
- 2701410 Battery charging unit
- 5805100 1 box of 100 Tonosafe prisms and 5 prism holders
- 2701421 Adaptor (UK)
- 2701422 Adaptor (EU)
- 2701423 Adaptor (US)

A1.1 Haag-Streit Goldmann measuring prisms

The original Haag-Streit measuring prisms are made of PMMA. Reusable measuring prisms must be cleaned and disinfected prior to every use. See the manufacturer's instruction for use.

A1.2 Tonosafe sterile disposable prisms

Tonosafe sterile disposable prisms have been developed to avoid the necessity of cleaning and disinfecting applanation tonometer prisms. Tonosafe is available in sterile packages with 100 disposable prisms.

Appendix B – Legal regulations

The Perkins tonometer was designed and developed to conform with EN 60601-1, EN ISO 8612 & ISO 15004-1 standards.

The Perkins tonometer satisfies the electromagnetic compatibility requirements of EN 60601-1-2. The device has been designed to maintain the emissions of electromagnetic interference at a level which does not exceed the statutory guidelines and which does not affect other devices in the vicinity. Statutory accident regulations are to be observed.

The LED's fulfil the requirements and have been tested to EN60825-1 Part 1: Safety of Laser Products: Equipment classification requirements under normal operating conditions and those of single fault failure. The tested product is a diagnostic device which includes LEDs as a method of illumination.

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© HAAG-STREIT UK 2026
Haag-Streit Deutschland GmbH
Haag-Streit AG